

SIRT6 Mouse mAb

Catalog No: #64238

Package Size: #64238-1 50ul #64238-2 100ul


Description

Product Name	SIRT6 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IF;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human SIRT6 (AA: 141-250) expressed in E. Coli.
Target Name	SIRT6
Other Names	SIR2L6
Accession No.	Q8N6T7
Calculated MW	39kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

ICC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Immunofluorescence analysis of HeLa cells using SIRT6 mouse mAb. Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin.



Immunofluorescence analysis of HeLa cells using SIRT6 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin. Secondary antibody from Fisher (Cat#: 35503)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.