

SIRT7 Mouse mAb

Catalog No: #64240

Package Size: #64240-1 50ul #64240-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	SIRT7 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human SIRT7 (AA: 1-105) expressed in E. Coli.
Target Name	SIRT7
Other Names	SIR2L7
Accession No.	Q9NRC8
Calculated MW	44.9kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC: 1/200 - 1/1000 ELISA: 1/10000


Images



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Immunohistochemical analysis of paraffin-embedded human liver cancer tissues using SIRT7 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human ovarian cancer tissues using SIRT7 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.