

SLC27A2 Mouse mAb

Catalog No: #64243

Package Size: #64243-1 50ul #64243-2 100ul


Description

Product Name	SLC27A2 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human SLC27A2 (AA: 346-405) expressed in E. Coli.
Target Name	SLC27A2
Other Names	VLCS; FATP2; VLACS; ACSVL1; FACVL1; hFACVL1; HsT17226
Accession No.	O14975
Calculated MW	70.3kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC: 1/200 - 1/1000 ELISA: 1/10000


Images



Black line: Control Antigen (100 ng); Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Immunohistochemical analysis of paraffin-embedded human liver cancer tissues using SLC27A2 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human esophageal cancer tissues using SLC27A2 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.