

SLC27A5 Mouse mAb

Catalog No: #64244

Package Size: #64244-1 50ul #64244-2 100ul


Description

Product Name	SLC27A5 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC
Species Reactivity	Human,Mouse
Immunogen Description	Purified recombinant fragment of human SLC27A5 (AA: 508-570) expressed in E. Coli.
Target Name	SLC27A5
Other Names	BAL; ACSB; BACS; FATP5; ACSVL6; FACVL3; FATP-5; VLACSR; VLCSH2; VLCS-H2
Accession No.	Q9Y2P5
Calculated MW	75.4kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200 - 1/1000ELISA:1/10000


Images




Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Western blot analysis using SLC27A5 mouse mAb against 3T3L1 (1) and NIH3T3 (2) cell lysate.



Immunohistochemical analysis of paraffin-embedded human liver cancer tissues using SLC27A5 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human stomach tissues using SLC27A5 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.