

SLC45A3 Mouse mAb

Catalog No: #64246

Package Size: #64246-1 50ul #64246-2 100ul


Description

Product Name	SLC45A3 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2b
Applications	WB;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human SLC45A3 (AA: 354-553) expressed in E. Coli.
Target Name	SLC45A3
Other Names	PRST; IPCA6; IPCA-2; IPCA-6; IPCA-8; PCANAP2; PCANAP6; PCANAP8
Accession No.	Q96JT2
Calculated MW	59 kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000FC:1/200 - 1/400ELISA:1/10000


Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using SLC45A3 mouse mAb against LNCap (1), PC-3 (2), HepG2 (3), SMMC-7721 (4), and DU145 (5) cell lysate.



Flow cytometric analysis of PC-3 cells using SLC45A3 mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.