

SMAD1 Mouse mAb

Catalog No: #64251

Package Size: #64251-1 50ul #64251-2 100ul


Description

Product Name	SMAD1 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Purified recombinant fragment of human SMAD1 (AA: 1-110) expressed in E. Coli.
Target Name	SMAD1
Other Names	BSP1; JV41; BSP-1; JV4-1; MADH1; MADR1
Accession No.	Q15797
Calculated MW	52.3kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000ELISA:1/10000

Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using SMAD1 mouse mAb against Hela (1), HepG2 (2), MCF-7 (3), C2C12 (4), A431 (5), A549 (6), LNCap (7), NIH/3T3 (8) and PC-12 (9) cell lysate.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.