

SMARCA2 Mouse mAb

Catalog No: #64254

Package Size: #64254-1 50ul #64254-2 100ul


Description

Product Name	SMARCA2 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC;IF;FC
Species Reactivity	Human
Target Name	SMARCA2
Other Names	BIS;BRM;SNF2;SWI2;hBRM;NCBRS;Sth1p;BAF190;SNF2L2;SNF2LA;hSNF2a
Accession No.	P51531
Calculated MW	181.3kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500-1/1000IHC:1/200-1/400ICC:1/25-1/50FC:1/200 - 1/400


Images




Western blot analysis using SMARCA2 mouse mAb against HeLa(1),HEK293(2),HepG2(3),A431(4)cell lysate.




Confocal Immunofluorescence analysis of NIH3T3 cells using SMARCA2 mouse mAb (green). Red: Actin filaments have been labeled with Alexa Fluor-555 phalloidin. Blue: DRAQ5 fluorescent DNA dye.



Flow cytometric analysis of HeLa cells using SMARCA2 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human ovarian cancer tissues using SMARCA2 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human rectal cancer tissues using SMARCA2 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.