

SNAI1 Mouse mAb

Catalog No: #64257

Package Size: #64257-1 50ul #64257-2 100ul


Description

Product Name	SNAI1 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human SNAI1 expressed in E. Coli.
Target Name	SNAI1
Other Names	SNA; SNAH; SLUGH2; dJ710H13.1; SNAI1
Accession No.	O95863
Calculated MW	29kDa
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500-1/1000ELISA:1/10000

Images



Red: Control Antigen (100ng); Purple: Antigen (10ng); Green: Antigen (50ng); Blue: Antigen (100ng);



Western blot analysis using SNAI1 mouse mAb against NTERA-2 cell lysate.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.