

SORL1 Mouse mAb

Catalog No: #64265

Package Size: #64265-1 50ul #64265-2 100ul


Description

Product Name	SORL1 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC;IF
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human SORL1 expressed in E. Coli.
Target Name	SORL1
Other Names	SORL1
Accession No.	Q92673
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200 - 1/1000ICC:1/200 - 1/1000ELISA:1/10000


Images



Confocal Immunofluorescence analysis of PANC-1 cells using SORL1 mouse mAb (green). Red: Actin filaments have been labeled with Alexa Fluor-555 phalloidin. Blue: DRAQ5 fluorescent DNA dye.



Immunohistochemical analysis of paraffin-embedded human cerebrum tissues using SORL1 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human brain, cortex tissues using SORL1 mouse mAb.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.