

SOX11 Mouse mAb

Catalog No: #64268

Package Size: #64268-1 50ul #64268-2 100ul


Description

Product Name	SOX11 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human SOX11 (AA: 1-250) expressed in E. Coli.
Target Name	SOX11
Other Names	CSS9; MRD27
Accession No.	P35716
Calculated MW	46.7kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200-1/1000FC:1/200-1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng)




Western blot analysis using SOX11 mouse mAb against Y-79 (1) cell lysate.



Flow cytometric analysis of Hela cells using SOX11 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human cervical cancer tissues using SOX11 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human bladder cancer tissues using SOX11 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.