

SP17 Mouse mAb

Catalog No: #64271

Package Size: #64271-1 50ul #64271-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	SP17 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human SP17 (AA: 1-152) expressed in E. Coli.
Target Name	SP17
Other Names	CT22; SPA17; SP17-1
Accession No.	Q15506
Calculated MW	17.4 kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200-1/1000FC:1/200-1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng)




Flow cytometric analysis of SK-OV-3 cells using SP17 mouse mAb (green) and negative control (red).



Flow cytometric analysis of Hela cells using SP17 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human rectal cancer tissues using SP17 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human testis tissues using SP17 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.