

## SPI1 Mouse mAb

Catalog No: #64273

Package Size: #64273-1 50ul #64273-2 100ul


## Description

|                       |                                                                                          |
|-----------------------|------------------------------------------------------------------------------------------|
| Product Name          | SPI1 Mouse mAb                                                                           |
| Host Species          | Mouse                                                                                    |
| Clonality             | Monoclonal                                                                               |
| Isotype               | Mouse IgG1                                                                               |
| Applications          | FC                                                                                       |
| Species Reactivity    | Human                                                                                    |
| Immunogen Description | Purified recombinant fragment of human SPI1 (AA: 124-271) expressed in E. Coli.          |
| Target Name           | SPI1                                                                                     |
| Other Names           | OF; PU.1; SFPI1; SPI-1; SPI-A                                                            |
| Accession No.         | P17947                                                                                   |
| Calculated MW         | 42kDa                                                                                    |
| Formulation           | Ascitic fluid containing 0.03% sodium azide.                                             |
| Storage               | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |


## Application Details

FC:1/200 - 1/400ELISA:1/10000

## Images



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Flow cytometric analysis of HeLa cells using SPI1 mouse mAb (green) and negative control (red).

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.