

SST Mouse mAb

Catalog No: #64277

Package Size: #64277-1 50ul #64277-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	SST Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC;FC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Purified recombinant fragment of human SST (AA: 1-116) expressed in E. Coli.
Target Name	SST
Other Names	SMST
Accession No.	P61278
Calculated MW	12.7kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);




Flow cytometric analysis of HepG2 cells using SST mouse mAb (green) and negative control (red).




Immunohistochemical analysis of paraffin-embedded human pancreas tissues using SST mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human lung cancer tissues using SST mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded Mouse kidney using SST mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded Rat kidney using SST mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.