

SSTR3 Mouse mAb

Catalog No: #64279

Package Size: #64279-1 50ul #64279-2 100ul


Description

Product Name	SSTR3 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;FC
Species Reactivity	Human,Rat
Immunogen Description	Purified recombinant fragment of human SSTR3 (AA: 1-43) expressed in E. Coli.
Target Name	SSTR3
Other Names	SS3R; SS3-R; SS-3-R; SSR-28
Accession No.	P32745
Calculated MW	45.8kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000FC:1/200 - 1/400ELISA:1/10000


Images



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Western blot analysis using SSTR3 mouse mAb against Hela (1), PANC-1 (2), PC-12 (3),SK-N-SH (4), U937 (5) and HepG2 (6) cell lysate.



Flow cytometric analysis of HeLa cells using SSTR3 mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.