

STAT5B Mouse mAb

Catalog No: #64282

Package Size: #64282-1 50ul #64282-2 100ul


Description

Product Name	STAT5B Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human STAT5B expressed in E. Coli.
Target Name	STAT5B
Other Names	STAT5; STAT5B
Accession No.	P51692
Calculated MW	90kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB:1/500 - 1/2000ELISA:1/10000

Images



Western blot analysis using STAT5B mouse mAb against Hela (1), K562 (2), NIH/3T3 (3), C6 (4), HEK293 (5), Jurkat (6) and HL-60 (7) cell lysate.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.