

STAT6 Mouse mAb

Catalog No: #64283

Package Size: #64283-1 50ul #64283-2 100ul


Description

Product Name	STAT6 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IF
Species Reactivity	Human,Mouse,Rat,Monkey
Immunogen Description	Purified recombinant fragment of human STAT6 expressed in E. Coli.
Target Name	STAT6
Other Names	STAT6B; STAT6C; D12S1644; IL-4-STAT; STAT6
Accession No.	P42226
Calculated MW	94kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000ICC:1/200 - 1/1000ELISA:1/10000

Images



Western blot analysis using STAT6 mouse mAb against HEK293 (1), NIH/3T3 (2), MCF-7 (3), Raw246.7 (4) and PC-12 (5) cell lysate.



Immunofluorescence analysis of HeLa cells using STAT6 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor-555 phalloidin

Note: This product is for in vitro research use only and is not intended for use in humans or animals.