

STK11 Mouse mAb

Catalog No: #64284

Package Size: #64284-1 50ul #64284-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	STK11 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;FC
Species Reactivity	Human,Mouse,Monkey
Immunogen Description	Purified recombinant fragment of human STK11 expressed in E. Coli.
Target Name	STK11
Other Names	PJS; LKB1; STK11
Accession No.	Q15831
Calculated MW	54kDa
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000FC:1/200 - 1/400ELISA:1/10000

Images



Western blot analysis using STK11 mouse mAb against NIH/3T3 (1),Raw246.7 (2), COS7 (3), Jurkat (4), HEK293 (5) and A431 (6) cell lysate.



Flow cytometric analysis of K562 cells using STK11 mouse mAb (green) and negative control (purple).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.