

SYK Mouse mAb

Catalog No: #64292

Package Size: #64292-1 50ul #64292-2 100ul


Description

Product Name	SYK Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human SYK (AA: 217-356) expressed in E. Coli.
Target Name	SYK
Other Names	p72-Syk
Accession No.	P43405
Calculated MW	72kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000FC:1/200 - 1/400ELISA:1/10000


Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using SYK mouse mAb against Ramos (1), HEK293 (2), Raji (3), and A431 (4) cell lysate.



Flow cytometric analysis of HeLa cells using SYK mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.