

SYP Mouse mAb

Catalog No: #64294

Package Size: #64294-1 50ul #64294-2 100ul


Description

Product Name	SYP Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC
Species Reactivity	Human,Rat
Immunogen Description	Purified recombinant fragment of SYP expressed in E. Coli.
Target Name	SYP
Other Names	MRXSYP; SYP
Accession No.	P08247
Calculated MW	33.8kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000ELISA:1/10000

Images



Western blot analysis using SYP mouse mAb against rat brain tissue lysate.



Immunohistochemical analysis of paraffin-embedded human Small Intestine, myenteric plexus tissues using anti-SYP mouse mAb

Note: This product is for in vitro research use only and is not intended for use in humans or animals.