

TCF3 Mouse mAb

Catalog No: #64301

Package Size: #64301-1 50ul #64301-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TCF3 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human TCF3 expressed in E. Coli.
Target Name	TCF3
Other Names	E2A; ITF1; VDIR; bHLHb21; MGC129647; MGC129648; TCF3
Accession No.	P15923
Calculated MW	67.6kDa
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000FC:1/200 - 1/400ELISA:1/10000

Images



Western blot analysis using TCF3 mouse mAb against A549 (1), A431 (2), Hela (3), PANC-1 (4) and PC-3 (5) cell lysate.



Flow cytometric analysis of A549 cells using anti-TCF3 mAb (green) and negative control (purple).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.