

TCL1A Mouse mAb

Catalog No: #64304

Package Size: #64304-1 50ul #64304-2 100ul


Description

Product Name	TCL1A Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human TCL1A (AA: 10-104) expressed in E. Coli.
Target Name	TCL1A
Other Names	TCL1
Accession No.	P56279
Calculated MW	13.5kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

FC:1/200 - 1/400ELISA:1/10000

Images



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Flow cytometric analysis of HeLa cells using TCL1A mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.