

TERT Mouse mAb

Catalog No: #64307

Package Size: #64307-1 50ul #64307-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TERT Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human TERT (AA: 1029-1132) expressed in E. Coli.
Target Name	TERT
Other Names	TP2; TRT; CMM9; EST2; TCS1; hTRT; DKCA2; DKCB4; hEST2; PFBMFT1
Accession No.	O14746
Calculated MW	127kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Flow cytometric analysis of HeLa cells using TERT mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human esophageal cancer tissues using TERT mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human rectum cancer tissues using TERT mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.