

TFAP2C Mouse mAb

Catalog No: #64309

Package Size: #64309-1 50ul #64309-2 100ul


Description

Product Name	TFAP2C Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human TFAP2C (AA: 341-450) expressed in E. Coli.
Target Name	TFAP2C
Other Names	ERF1; TFAP2G; hAP-2g; AP2-GAMMA
Accession No.	Q92754
Calculated MW	49.2kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000ELISA:1/10000


Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using TFAP2C mouse mAb against SK-Br-3 (1) cell lysate.



Immunohistochemical analysis of paraffin-embedded human esophageal cancer tissues using TFAP2C mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.