

TFRC Mouse mAb

Catalog No: #64311

Package Size: #64311-1 50ul #64311-2 100ul

Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com


Description

| | |
|-----------------------|--|
| Product Name | TFRC Mouse mAb |
| Host Species | Mouse |
| Clonality | Monoclonal |
| Isotype | Mouse IgG1 |
| Applications | WB |
| Species Reactivity | Human,Mouse,Monkey,Rat |
| Immunogen Description | Purified recombinant fragment of human TFRC (AA: 608-727) expressed in E. Coli. |
| Target Name | TFRC |
| Other Names | T9; TR; TFR; p90; CD71; TFR1; TRFR |
| Accession No. | P02786 |
| Calculated MW | 85kDa |
| Formulation | Purified antibody in PBS with 0.05% sodium azide |
| Storage | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |


Application Details

WB:1/500 - 1/2000ELISA:1/10000

Images



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Western blot analysis using TFRC mouse mAb against Jurkat (1), Hela (2), K562 (3), Cos7 (4), MCF-7 (5), PC-12 (6), NIH/3T3 (7), HEK293 (8), RAJI (9) cell lysate.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.