

THBS1 Mouse mAb

Catalog No: #64312

Package Size: #64312-1 50ul #64312-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	THBS1 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human THBS1 (AA: 750-850) expressed in E. Coli.
Target Name	THBS1
Other Names	TSP; THBS; TSP1; TSP-1; THBS-1
Accession No.	P07996
Calculated MW	129kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

FC:1/200 - 1/400ELISA:1/10000

Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of Hela cells using THBS1 mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.