

Thioredoxin (TRX) Mouse mAb

Catalog No: #64313

Package Size: #64313-1 50ul #64313-2 100ul


Description

Product Name	Thioredoxin (TRX) Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB
Species Reactivity	Tag
Immunogen Description	Purified recombinant fusion protein with Thioredoxin (TRX) tag.
Target Name	Thioredoxin (TRX)
Other Names	Thioredoxin
Calculated MW	12kDa
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB:1/500 - 1/2000IHC:1/100 - 1/1000ICC:1/100-1/200FC:1/200 - 1/400ELISA:1/10000

Images



Western blot analysis using Trx mouse mAb against various fusion protein with Trx tag.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.