

TIMD4 Mouse mAb

Catalog No: #64320

Package Size: #64320-1 50ul #64320-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TIMD4 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human TIMD4 (AA: extra 25-314) expressed in E. Coli.
Target Name	TIMD4
Other Names	TIM4; SMUCKLER
Accession No.	Q96H15
Calculated MW	41.6kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000FC:1/200 - 1/400ELISA:1/10000


Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using TIMD4 mouse mAb against Jurkat (1), HEK293 (2), HL-60 (3), and Hela (4) cell lysate.



Flow cytometric analysis of Jurkat cells using TlMD4 mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.