

Tlr2 Mouse mAb

Catalog No: #64323

Package Size: #64323-1 50ul #64323-2 100ul


Description

Product Name	Tlr2 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of mouse Tlr2 (AA: 628-777) expressed in E. Coli.
Target Name	Tlr2
Other Names	Ly105
Accession No.	Q9QUN7
Calculated MW	90kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000ELISA:1/10000

Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using Tlr2 mouse mAb against Ramos (1), U937 (2), Jurkat (3), and Hela (4) cell lysate.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.