

TNFRSF10D Mouse mAb

Catalog No: #64330

Package Size: #64330-1 50ul #64330-2 100ul


Description

Product Name	TNFRSF10D Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB
Species Reactivity	Human,Monkey
Immunogen Description	Purified recombinant fragment of human TNFRSF10D (AA: extra 56-211) expressed in E. Coli.
Target Name	TNFRSF10D
Other Names	DCR2; CD264; TRUNDD; TRAILR4; TRAIL-R4
Accession No.	Q9UBN6
Calculated MW	41.8kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000ELISA:1/10000

Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using TNFRSF10D mouse mAb against Hela (1), L-02 (2), A549 (3), HepG2 (4), and COS7 (5) cell lysate.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.