

TNFRSF11A Mouse mAb

Catalog No: #64332

Package Size: #64332-1 50ul #64332-2 100ul



Signalway Antibody

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	TNFRSF11A Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;FC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Purified recombinant fragment of human TNFRSF11A (AA: extra 30-212) expressed in E. Coli.
Target Name	TNFRSF11A
Other Names	FEO; OFE; ODFR; OSTs; PDB2; RANK; CD265; OPTB7; TRANCER; LOH18CR1; TRANCE-R
Accession No.	Q9Y6Q6
Calculated MW	66kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000FC:1/200 - 1/400ELISA:1/10000


Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using TNFRSF11A mouse mAb against Mouse kidney (1) and Rat kidney (2) cell lysate.



Flow cytometric analysis of HL-60 cells using TNFRSF11A mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.