

TNFRSF11B Mouse mAb

Catalog No: #64335

Package Size: #64335-1 50ul #64335-2 100ul


Description

Product Name	TNFRSF11B Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human TNFRSF11B expressed in E. Coli.
Target Name	TNFRSF11B
Other Names	OPG; TR1; OCIF; MGC29565; TNFRSF11B
Accession No.	O00300
Calculated MW	60kDa
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

IHC: 1/200 - 1/1000 ELISA: 1/10000

Images



Immunohistochemical analysis of paraffin-embedded human mammary cancer tissues (left) and lung cancer tissues (right) using TNFRSF11B mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.