

TNFRSF25 Mouse mAb

Catalog No: #64339

Package Size: #64339-1 50ul #64339-2 100ul


Description

Product Name	TNFRSF25 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human TNFRSF25 (AA: extra(25-199)) expressed in E. Coli.
Target Name	TNFRSF25
Other Names	DR3; TR3; DDR3; LARD; APO-3; TRAMP; WSL-1; GEF720; WSL-LR; PLEKHG5; TNFRSF12
Accession No.	Q93038
Calculated MW	45.4kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200-1/1000FC:1/200-1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng)



Flow cytometric analysis of THP-1 cells using TNFRSF25 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human cervical cancer tissues using TNFRSF25 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human stomach cancer tissues using TNFRSF25 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.