

TNNT2 Mouse mAb

Catalog No: #64344

Package Size: #64344-1 50ul #64344-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TNNT2 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC;FC
Species Reactivity	Human,Mouse
Immunogen Description	Purified recombinant fragment of human TNNT2 (AA: full 1-295) expressed in E. Coli.
Target Name	TNNT2
Other Names	CMH2; RCM3; TnTC; cTnT; CMD1D; CMPD2; LVNC6
Accession No.	P45379
Calculated MW	35.9kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using TNNT2 mouse mAb against mouse heart (1) cell lysate.



Flow cytometric analysis of Hale cells using TNNT2 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human myocardium tissues using TNNT2 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.