

TOP1 Mouse mAb

Catalog No: #64345

Package Size: #64345-1 50ul #64345-2 100ul


Description

Product Name	TOP1 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IF;FC
Species Reactivity	Human,Mouse
Immunogen Description	Purified recombinant fragment of human TOP1 expressed in E. Coli.
Target Name	TOP1
Other Names	TOPI
Accession No.	P11387
Calculated MW	90.7kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000 ICC:1/100 - 1/400 FC:1/200 - 1/400


Images



Western blot analysis using TOP1 mouse mAb against HUVEC (1), Jukat (2), MCF-7 (3), PC-12 (4), Hela (5), NIH/3T3 (6) and HCT116 (7) cell lysate.



Immunofluorescence analysis of Hela cells using TOP1 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor-555 phalloidin. Secondary antibody from Fisher (Cat#: 35503)



Flow cytometric analysis of HeLa cells using TOP1 mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.