

TOP2A Mouse mAb

Catalog No: #64346

Package Size: #64346-1 50ul #64346-2 100ul


Description

Product Name	TOP2A Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human TOP2A (AA: 1100-1530) expressed in E. Coli.
Target Name	TOP2A
Other Names	TOP2; TP2A
Accession No.	P11388
Calculated MW	174.4kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000ELISA:1/10000


Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using TOP2A mouse mAb against HepG2 (1), Hela (2), Jurkat (3), and K562 (4) cell lysate.



Immunohistochemical analysis of paraffin-embedded human rectum cancer tissues using TOP2A mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.