

TP63 Mouse mAb

Catalog No: #64349

Package Size: #64349-1 50ul #64349-2 100ul


Description

Product Name	TP63 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human TP63 (AA: 5-17) expressed in E. Coli.
Target Name	TP63
Other Names	TP63,AIS; KET; LMS; NBP; RHS; p40; p51; p63; EEC3; OFC8; p73H; p73L; SHFM4; TP53L; TP73L; p53CP; TP53CP; B(p51A); B(p51B)
Accession No.	Q9H3D4
Calculated MW	76.8kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200-1/1000FC:1/200-1/400ELISA:1/10000


Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng)



Flow cytometric analysis of Hela cells using TP63 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human bladder Cancer tissues using TP63 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.