

TRAF2 Mouse mAb

Catalog No: #64354

Package Size: #64354-1 50ul #64354-2 100ul


Description

Product Name	TRAF2 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2b
Applications	WB;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human TRAF2 (AA: 39-188) expressed in E. Coli.
Target Name	TRAF2
Other Names	TRAP; TRAP3; MGC:45012
Accession No.	Q12933
Calculated MW	55.8kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000FC:1/200 - 1/400ELISA:1/10000


Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using TRAF2 mouse mAb against MCF-7 (1), A431 (2), Hela (3), Jurkat (4), HEK293 (5), and Ramos (6) cell lysate.



Flow cytometric analysis of HeLa cells using TRAF2 mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.