

TRAFD1 Mouse mAb

Catalog No: #64355

Package Size: #64355-1 50ul #64355-2 100ul


Description

|                       |                                                                                          |
|-----------------------|------------------------------------------------------------------------------------------|
| Product Name          | TRAFD1 Mouse mAb                                                                         |
| Host Species          | Mouse                                                                                    |
| Clonality             | Monoclonal                                                                               |
| Isotype               | Mouse IgG1                                                                               |
| Applications          | WB;IF;FC                                                                                 |
| Species Reactivity    | Human                                                                                    |
| Immunogen Description | Purified recombinant fragment of human TRAFD1 (AA: 401-582) expressed in E. Coli.        |
| Target Name           | TRAFD1                                                                                   |
| Other Names           | FLN29                                                                                    |
| Accession No.         | O14545                                                                                   |
| Calculated MW         | 64.8kDa                                                                                  |
| Formulation           | Purified antibody in PBS with 0.05% sodium azide                                         |
| Storage               | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |


Application Details

ICC:1/100FC:1/200 - 1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Western blot analysis using TRAFD1 mouse mAb against HEK293 (1), Raji (2), and Jurkat (3) cell lysate.



Immunofluorescence analysis of HepG2 cells using TRAFD1 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye.



Flow cytometric analysis of HeLa cells using TRAFD1 mouse mAb (green) and negative control (purple).

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.