

Trim5a Mouse mAb

Catalog No: #64361

Package Size: #64361-1 50ul #64361-2 100ul


Description

Product Name	Trim5a Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human trim5 alpha expressed in E. Coli.
Target Name	Trim5a
Other Names	RNF88; TRIM5alpha
Accession No.	Q9C035
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000ELISA:1/10000

Images



Western blot analysis using Trim5a mouse mAb against human breast carcinoma tissue lysate.



Note: This product is for in vitro research use only and is not intended for use in humans or animals.