

TRIP6 Mouse mAb

Catalog No: #64362

Package Size: #64362-1 50ul #64362-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TRIP6 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human TRIP6 (AA: 107-291) expressed in E. Coli.
Target Name	TRIP6
Other Names	OIP1; OIP-1; ZRP-1; TRIP-6; TRIP6i2
Accession No.	Q15654
Calculated MW	50.3kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000FC:1/200 - 1/400ELISA:1/10000


Images



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Western blot analysis using TRIP6 mouse mAb against K562 and A431 (2) cell lysate.



Flow cytometric analysis of K562 cells using TRIP6 mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.