

troponin T2 Mouse mAb

Catalog No: #64363

Package Size: #64363-1 50ul #64363-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	troponin T2 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human troponin T2 expressed in E. Coli.
Target Name	troponin T2
Other Names	CMH2; RCM3; TnTC; cTnT; CMPD2; LVNC6
Accession No.	P45379
Calculated MW	36kDa
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

FC:1/200 - 1/400ELISA:1/10000

Images



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Flow cytometric analysis of MCF-7 cells using troponin T2 mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.