

TSLPR Mouse mAb

Catalog No: #64366

Package Size: #64366-1 50ul #64366-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TSLPR Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human TSLPR (AA: extra(23-231)) expressed in E. Coli.
Target Name	TSLPR
Other Names	CRL2; TSLPR; CRLF2Y
Accession No.	Q9HC73
Calculated MW	42kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200-1/1000FC:1/200-1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng)



Flow cytometric analysis of THP-1 cells using TSLPR mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human bladder cancer tissues using TSLPR mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human rectal cancer tissues using TSLPR mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.