

TSPAN33 Mouse mAb

Catalog No: #64368

Package Size: #64368-1 50ul #64368-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TSPAN33 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2b
Applications	IHC;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human TSPAN33 (AA:Extra(46-64)+Extra(118-235) expressed in E. Coli.
Target Name	TSPAN33
Other Names	PEN; PEN.; TSPAN-33
Accession No.	Q86UF1
Calculated MW	32kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of Jurkat cells using TSPAN33 mouse mAb (green) and negative control (red).



Flow cytometric analysis of MOLT4 cells using TSPAN33 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human endometrial cancer tissues using TSPAN33 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.