

Tyro3 Mouse mAb

Catalog No: #64375

Package Size: #64375-1 50ul #64375-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Tyro3 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IF
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of Tyro3 (aa138-321) expressed in E. Coli.
Target Name	Tyro3
Other Names	BYK; Brt; Dtk; RSE; Sky; Tif; FLJ16467
Accession No.	Q06418
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

ICC:1/200 - 1/1000ELISA:1/10000

Images



Immunofluorescence analysis of Hela (Left) and MCF-7 (Right) cells using Tyro3 mouse mAb (green).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.