

TYRO3 Mouse mAb

Catalog No: #64376

Package Size: #64376-1 50ul #64376-2 100ul


Description

Product Name	TYRO3 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB
Species Reactivity	Human
Immunogen Description	Purified recombinant extracellular fragment of human TYRO3 fused with hIgGFc tag expressed in HEK293 cell line.
Target Name	TYRO3
Other Names	BYK; Brt; Dtk; RSE; Sky; Tif
Accession No.	Q06418
Calculated MW	96.9kDa
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

ELISA:1/10000

Images



Western blot analysis using TYRO3 mouse mAb against extracellular domain of human TYRO3 (aa41-429).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.