

TYRO3 Mouse mAb

Catalog No: #64377

Package Size: #64377-1 50ul #64377-2 100ul


Description

Product Name	TYRO3 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human TYRO3 (AA: extra 230-429) expressed in mammalian.
Target Name	TYRO3
Other Names	BYK; Dtk; RSE; Rek; Sky; Tif; Etk-2
Accession No.	Q06418
Calculated MW	96.9kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of U937 cells using TYRO3 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human liver cancer tissues using TYRO3 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human brain tissues using TYRO3 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.