

UBB Mouse mAb

Catalog No: #64378

Package Size: #64378-1 50ul #64378-2 100ul


Description

Product Name	UBB Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human UBB expressed in E. Coli.
Target Name	UBB
Other Names	UBC; UBA52; RPS27A
Accession No.	P0CG47
Calculated MW	26kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000FC:1/200 - 1/400ELISA:1/10000


Images



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Western blot analysis using UBB mouse mAb against NIH/3T3 (1) and HeLa (2) cell lysate.



Flow cytometric analysis of Hela cells using UBB mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.