

UBE2C Mouse mAb

Catalog No: #64379

Package Size: #64379-1 50ul #64379-2 100ul


Description

Product Name	UBE2C Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2a
Applications	WB;IHC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human UBE2C (AA: FULL(1-179)) expressed in E. Coli.
Target Name	UBE2C
Other Names	UBCH10; dJ447F3.2
Accession No.	O00762
Calculated MW	19.7kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200 - 1/1000ELISA:1/10000


Images




Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Western blot analysis using UBE2C mouse mAb against Hela (1), and Raji (2) cell lysate.



Immunohistochemical analysis of paraffin-embedded human liver cancer tissues using UBE2C mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human esophageal cancer tissues using UBE2C mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.