

UCP3 Mouse mAb

Catalog No: #64382

Package Size: #64382-1 50ul #64382-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	UCP3 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2a
Applications	IF
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human UCP3 (AA: 1-113 and 217-312) expressed in E. Coli.
Target Name	UCP3
Other Names	SLC25A9
Accession No.	P55916
Calculated MW	34.2kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

ICC:1/200 - 1/1000ELISA:1/10000

Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Immunofluorescence analysis of HL-7702 cells using UCP3 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin. Secondary antibody from Fisher (Cat#: 35503)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.