

UL37 Mouse mAb

Catalog No: #64384

Package Size: #64384-1 50ul #64384-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	UL37 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC;IF;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of Human herpesvirus UL37 (AA: 970-1119) expressed in E. Coli.
Target Name	UL37
Accession No.	Q96RL1
Calculated MW	120.5kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200 - 1/1000 ICC:1/50 - 1/250 FC:1/200 - 1/400 ELISA:1/10000


Images




Black line: Control Antigen (100 ng); Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of MCF-7 cells using UL37 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human bladder cancer tissues using UL37 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human esophageal cancer tissues using UL37 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.